

Quick Start Guide

Wireless Pager Standalone Kit

Applies to product codes:

R14-0001 | R14-0002 | R14-0003 | R14-0004 | R14-0005 |
R14-0006 | R14-0008

| System Contents | |
|-----------------|--|
| Monitor | RG41 Wireless Jill |
| Sensor | RG49 Wireless Seatgard RG47 Wireless Bedgard RG38 Wireless Floorgard / RG59 Wireless Floorgard RG39 Wireless Under-mattress Pad |
| Transmitter | RG9 Wireless Pad Transmitter |

Sensor Placement

RG47 Bedgard: Place on top of the mattress and under the bed sheet. The pad can be placed anywhere between the shoulders to the hips depending on the condition of the patient and how early you want the alert.

RG39 Under-mattress: Place under the mattress across the width of the bed at the lower back level. On profiling beds make sure it will not entrap the mechanism when the bed is being adjusted.

For further guidance on placement, please refer to the pad positioning flyer.

RG49 Seatgard: Place mat on the relevant seat on top of or under a cushion, giving consideration to the position of the monitor and the least intrusive method of installing the equipment.

RG38/RG59 Floorgard: Place mat at or near to the area to be monitored such as a doorway or side of the bed.

Receiver/Monitor Placement

The Jill pager can be wall mounted using the supplied bracket, worn around the neck on a lanyard or clipped onto a pocket.

General Operation

The pager device can receive signals from a sensor and reset alarms by pressing and holding the cancel button. Multiple pagers can also be used for this set-up.



Quick Start Guide

Wireless Pager Standalone Kit

Follow the steps below to set up your equipment:

- 1.** Pull out the battery tab to power up the Wireless Jill or insert 2 AA batteries.
- 2.** Switch the pager onto either **VIBRATE & TONE** or **VIBRATE ONLY**.
- 3.** Ensure that the sensor pad is placed on a smooth flat surface and the pad is not creased or folded in any way.
- 4.** Check that you have positioned your sensor pad correctly and switch the transmitter on (inside pouch).
- 5.** When you first switch on the pad transmitter, apply pressure to the sensor pad for at least 15 secs to wake up the transmitter from sleep mode. Failure to do this will cause the Companion to loose connection from the transmitter and trigger an audible alert. This only applies when you switch a transmitter from OFF to ON for the first time.
- 6.** To test, activate the sensor to confirm that it is alarming and resetting correctly.

When you are satisfied with the points above please complete the daily checks sheet provided before putting the system into use.

**Ramblegard Technical Support Helpline:
028 9448 3979**

