

# Companion Monitor

## Re-programming Instructions

Follow the steps below to re-program a sensor pad into the Companion monitor:

- 1.** To re-program a sensor pad into the Companion, go to the Companion and hold the **PROGRAM** button (on the back) and **RESET** button together until it beeps and '**01**' is displayed on the screen.
- 2.** If there is currently a pad paired into '**01**', the pad status light will be illuminated.
- 3.** If illuminated, hold **RESET** until the monitor beeps and the pad status light goes out.
- 4.** Your Companion is now ready to pair with a pad. Go to the transmitter on the pad and switch from **OFF** to **ON**. The Companion will beep, confirming it has received the signal and the pad status will be illuminated.
- 5.** Should you have a second pad to program, press **RESET** once and '**02**' is displayed on the screen. Repeat steps 2-4.
- 6.** To take out of program mode, press the **PROGRAM** button on the back of the monitor until it beeps. *Your system is now set up and ready for use.*

### PLEASE NOTE:

1. We recommend always clearly marking a matching set with corresponding labels or numbers so that you keep the paired pad and pager together.
2. There is a tracking option on the Companion, if ON, it constantly looks for the pad it is paired to, which can cause confusion. We recommend switching this OFF unless using the system with a pendant.

### Ramblegard Technical Support Helpline:

**+44 (0)28 9448 3979**

