

# Corded Ramblegard Monitor



Picture for illustrative purposes only

## Product Code:

RG1

## Description:

The Ramblegard Corded Monitor will work in conjunction with corded bed/seat or floor mat and alert accordingly.

## Features:

- Visual alert LED activates during alarm
- Reset Button
- Choice of tone and volume control for alarm
- Port to connect with nurse call system (NO)
- Operates on 3 x AA batteries
- Low battery and nurse call plug removal alert
- Pad lost alert
- Delay Feature - 2sec, 5min, 15min.

## Works in conjunction with:

- RG12 - Wired Bedgard
- RG10 - Wired Seatgard
- RG15 - Wired Floorgard

## Technical Specifications:

Power Supply :	3 x AA-type alkaline batteries
Battery Life :	6 months (Nominal)
Dimensions :	(LxWxD): 120 x 77 x 34mm
Housing material :	ABS (UL-94HB, flammability grade)
Net Weight:	150g Operating Temp. 5°C ~ 40°C
Rel. Humidity:	< 85%

**It stands out from its competitors,  
it is effective and it's economical.**

**Ramblegard**<sup>®</sup>  
stay a step ahead

[www.ramblegard.com](http://www.ramblegard.com)